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EN

**Instruction for Use****RABBIT COAGULASE PLASMA (LYOPHILIZED)****Intended use**

The rabbit coagulase plasma (lyophilized) is an *in vitro* diagnostic medical device used to identify coagulase-positive staphylococci (including *Staphylococcus aureus*). It is a qualitative test used to determine free and bound coagulase (clumping factor) produced by *Staphylococcus aureus*.

**Description**

Product for professional laboratory application only that should be used by qualified personnel following basic aseptic and antiseptic principles.

The plasma is in a glass vial with a stopper and a flip-off cap. One vial contains lyophilized rabbit plasma in an amount:

- 2 ml for POK02

- 5ml for POK05

The product must be rehydrated before use according to the instructions.

**Composition**

The plasma from rabbit blood collected with an CPDA or sodium citrate anticoagulant.

**Properties**

Coagulase is an enzyme produced by some strains of staphylococci - so-called coagulase-positive staphylococci (including *Staphylococcus aureus*). Coagulase in staphylococci occurs in two forms - free (released into the environment) and bound with the bacterial cell wall.

Free coagulase is detected in the tube test method by mixing rehydrated plasma in a tube with a 24-hour bacterial culture. The result is checked after the appropriate incubation time. The appearance of a clot in the tube indicates a positive result. In the tube test, coagulase protein is released into the suspension by coagulase-positive bacteria. Coagulase converts the fibrinogen present in the plasma into fibrin, which is then observed in the form of a clot.

On the other hand, bound coagulase (also known as clumping factor) is determined by the slide test. The slide test is performed in the following manner: a fresh bacterial culture taken from a non-selective agar medium is suspended in a drop of saline on a glass slide, and then with the use of a sterile loop rehydrated plasma is added to the mixture. The appearance of characteristic flocs after adding the rabbit plasma drop indicates a positive result. Clumping factor is a fibrinogen receptor on the surface of coagulase-positive bacterial cells. Fibrinogen, by converting to fibrin, forms threads between cells linking them in aggregates visible to the naked eye.

**Method****Preparation of samples**

Samples should be collected to sterile containers, transported to the laboratory and processed according to guidelines for the specific type of sample. Clearly defined, typical colonies should be selected. Uncontaminated bacterial cultures should be used for testing.

**Preparation of bacterial culture**

To test for free coagulase using the tube method, the selected strain should be inoculated to a liquid, growth-promoting medium such as enriched broth and incubated for 24 ± 2 hours at 37°C.

To test for bound coagulase using the slide method, the selected strain should be streaked on the appropriate solid (containing agar), growth-promoting medium and incubated for 24 hours at 37°C.

**Plasma preparation**

Before use, the plasma should be rehydrated by aseptically adding the volume of sterile deionized water indicated on the label (2 mL or 5 mL) to the vial. After adding the water, immediately begin gentle mixing and continue until the plasma is completely dissolved; do not leave the vial without mixing during rehydration.

**Plasma freezing:** If plasma freezing is planned, it is recommended to aliquot the material prior to freezing into small volumes, for example 0.5 mL (single-reaction volume). Frozen material should be thawed at 18–25°C. During thawing, gently invert the vial periodically to ensure uniform thawing. Do not leave the vial standing without mixing until completely thawed. To accelerate thawing, the closed vial may be gently warmed in the hand while mixing. After complete thawing, mix gently until a homogeneous solution is obtained.

In order to maintain the aseptic nature of the product (lack of which can prevent the test from being performed or the results from being properly interpreted), the preparation of the plasma and the performance of the test should take place in conditions that prevent contamination or pollution of the product.

**Procedure for determining free coagulase using the tube method**

1. Add 500 µl of rehydrated plasma to a sterile test tube. Add 500 µl of the strain from the broth culture.
2. Incubate test tubes for 24 ± 2 hours at 37 ± 1 °C.


**Reading and interpreting the results**

The results should be checked after 2, 4, 6 and 24 ± 2 hours from the start of incubation.

The appearance of a clot in the tube indicates a positive test result. The absence of a clot in the tube indicates a negative result.

! Do not shake or mix the contents of the tube (this may lead to the obtaining of false-negative results due to the mechanical breakdown of the clot).

Interpretation of the test results should be made according to the following scale:

- negative: no clot
- 1+ positive: weak uncompact clot
- 2+ positive: weak compact clot
- 3+ positive: strong compact clot
- 4+ positive: the clot occupies the whole volume of the tube and does not move when the tube is inverted.

**NOTE:** Some strains of *Staphylococcus aureus* produce fibrinolysin - an enzyme that breaks down clots. This may cause false-negative results after 24 ± 2 hours, therefore the tubes should be observed not only after 24 ± 2 hours, but also after 2, 4, 6 hours from the start of incubation.

*Staphylococcus aureus* strains cause coagulation of rabbit plasma within half an hour to 24 hours. Usually the plasma is completely clotted (solid form in the tube) but any degree of clotting should be considered as a positive result.

An instructional video with detailed information on how to perform the tube method test and how to interpret the obtained results is available on our website: <http://immunolab.com.pl/filmy>.

#### Procedure for determining bound coagulase (clumping factor) using glass slide method

1. Place two drops of physiological saline solution next to each other on a clean, dry glass slide.
2. In both drops suspend the same bacterial culture taken from a non-selective agar medium. The suspensions should be dense and homogeneous.
3. Add 10 µl of rabbit plasma to one of the drops with a sterile bacteriological loop and mix gently.
4. Gently rock the slide in a circular motion. Observe the result of the reaction.

#### Interpretation of the test results

If characteristic clumps (precipitates) appear within **60 seconds** of adding rabbit plasma, then the result should be considered positive.

**NOTE:** The negative control (bacterial suspension without the addition of rabbit plasma) should remain homogeneous. Otherwise, the tested strain shows autoagglutination, which means that the clumping factor cannot be determined. The reaction is more visible if it is observed over a dark background.

An instructional video with detailed information on how to perform the glass slide method test and how to interpret the obtained results is available on our website: <http://immunolab.com.pl/filmy>.

#### Quality control performed by user

During the application of plasma it is recommended to perform control tests using standard coagulase-positive and coagulase-negative strains in order to verify the correct performance of the method and the effectiveness of plasma, proceeding similarly as in the case of tested strains.

#### Limitations of the method

Besides *Staphylococcus aureus*, coagulase-positive bacteria that could react with rabbit plasma include other bacteria of the *Staphylococcus* genus, such as *S. anaerobius*, *S. delphini*, *S. hyicus*, *S. intermedius*, *S. lutrae* and *S. schleiferi* subsp. *coagulans*. However, these bacteria are less common than *Staphylococcus aureus* and are mainly responsible for infections in animals.

It should be noted that **not all** strains of *S. aureus* produce the coagulase enzyme and those are considered atypical. This means that some strains (especially clinical isolates or those showing methicillin resistance) may give a negative result in the test for the presence of free or bound coagulase.

For analysis, only pure bacterial cultures of staphylococci should be used. The presence of citrate may cause false positive reactions with other microorganisms (e.g. *Enterococcus faecalis*, *Serratia marcescens*, bacteria belonging to the *Pseudomonas* and *Streptococcus* genre).

#### Materials required to perform the assay not included in the kit

Dark background, physiological saline, sterile deionized water, droppers, test tubes, Petri dishes, glass slides, sterile pipette tips, container for contaminated materials, an instructional video (available on our website: <http://immunolab.com.pl/filmy>).

#### Validation and verification of the method

<b>sensitivity</b>	<b>free coagulase: 100% (38/38)</b> <b>clumping factor: 100% (35/35)</b> (number of true positive results/ sum of true positive and false negative results)
<b>specificity</b>	<b>free coagulase: 100% (14/14)</b> <b>clumping factor: 100% (16/16)</b> (number of true negative results/ sum of true negative and false negative results)
<b>repeatability</b>	<b>free coagulase: 100% (156/156)</b> <b>clumping factor: 100% (147/147)</b> (number of repeated compatible results/total number of results)

<b>reproducibility</b>	<b>free coagulase: 100% (149/149)</b> <b>clumping factor: 100% (141/141)</b> (number of reproduced compatible results/ total number of results)
------------------------	---

#### Storage conditions

- The **lyophilized plasma** should be stored at a temperature of **+2 to +8°C**.
- Transport for up to 3 days at temperatures not exceeding +15 to +25°C does not change the properties of the product.
- Do not use after the expiration date on the packaging.
- The expiration date applies to the product in intact packaging, stored in accordance with the instructions.
- The plasma is a reusable product.
- **Reconstituted plasma** should be stored at a temperature of **+2 to +8°C** no longer than **14 days** from the date of rehydration or at a temperature of **-20°C** no longer than **30 days**.
- **Repeated freezing and thawing is not permitted.**
- Protect from light.
- Do not use plasma if it is cloudy or has clotted. If the plasma is difficult to dissolve, do not use it.
- Do not use plasma with altered physico-chemical properties indicating microbiological contamination. Do not sterilize plasma prior to use.

#### Precautions

- The product does not pose a threat to the environment or human health if used in accordance with the instructions.
- Aseptic working conditions must be observed and precautions against microbiological hazards should be taken during all procedures.
- The instructions for use must be strictly followed.
- Any serious incident involving the device should be reported to the manufacturer and the competent authority of the Member State in which the user resides or has his place of business.

#### Precautions when working with infectious material:

- Rabbit plasma contains animal-derived material so should be treated as potentially infectious material and handled accordingly.
- The test may involve working with pathogenic bacteria so all necessary precautions for working with infectious microbiological material should be observed. Wear protective clothing and gloves. Collected material for testing and material on slides and in tubes should be protected from access by third parties due to the possibility of infection. Used slides, test tubes and other contaminated materials, should be sterilized in an autoclave or disposed of in accordance with the laboratory's standard procedures.

#### Contents of a single package

catalogue number	contents	number of tests- free coagulase/clumping factor
POK01	1x2ml	4 / 200
POK10	10x2ml	40 / 2000
POK05	1x5ml	10 / 500

#### Manufacturer

Research and Development Department of Salmonella Center "IMMUNOLAB" Ltd.  
Address: ul. Kładki 24, 80-822 Gdańsk, Poland  
Phone/Fax: 0048 58 781-44-91  
E-mail: [info@immunolab.com.pl](mailto:info@immunolab.com.pl)  
[www.immunolab.com.pl](http://www.immunolab.com.pl)

#### Explanation of the symbols used

Manufacturer

Expiration date : YYYY-XX-Z (year-month-day)

Batch code (Lot)

Consult instructions for use

Catalogue number

Storage temperature limit

In vitro diagnostic medical device

Transport temperature limit up to 3 days

Changes introduced in the current edition are highlighted in gray.